

HEALEON DUET™

PLATELET-RICH PLASMA (PRP)
PREPARATION SYRINGE



Designed for use within the practitioner's office, the **Healeon DUET** system with **Dial-In** technology delivers a tailored concentration approach for a variety of patient treatment applications.

REGULATORY STATUS:

DEVICE NAME: HEALEON DUET
CE MARK

Indications for Use: Healeon Duet is indicated to be used for the safe and rapid preparation of autologous platelet rich plasma (PRP) from a small sample of peripheral blood at the patient's point of care.

Learn more about the
Healeon DUET at
www.healeonmedical.com



HEALEON DUET™

PLATELET-RICH PLASMA (PRP)
PREPARATION SYRINGE



BIMINI HEALTH TECH
420 Stevens Ave., #220
Solana Beach, CA 92075 USA

+1-858-348-8050
info@biminihealthtech.com
www.biminihealthtech.com

Not for U.S.A. dissemination | 500-063-01 REV A.


Create Customized Platelet Concentration
Levels In Controlled Procedure-Specific Doses
With A **SINGLE PLATFORM**

ENGINEERED FOR PERFECTION

TAILORED FOR EVERY APPROACH




The **Healeon DUET** provides a customized system that allows the desired plasma and hematocrit fraction to be collected in a closed, simple and precise way, directly into the application syringe.



Dial-In™ technology for a tailored approach. You select the **OPTIMUM** PRP dose for your patients, whether in office or surgical settings, in a safe, closed-system device.

Operating similar to a standard syringe provides intuitive handling, with simple user training, allowing reproducibility of PRP quality that is tailored for any treatment approach.



DEVELOPED TO CARE FOR AND PROTECT THE PLATELETS

HEALEON DUET™

The **Healeon DUET** enables safe, quick, and easy separation of blood components into customizable platelet concentrations for specific procedures in a single platform. The most **VERSATILE** PRP on the market.



FULLY CUSTOMIZABLE PLASMA

- PRP (12 ml) or TOTAL PLASMA (from 22.5 ml Whole Blood)
- HIGH RECOVERY (mean 94.3% +/- 4.8% SD)
- TUNABLE HEMATOCRIT LEVELS



HIGHEST OUTPUT/PURITY

- HIGH PURITY (>98% of erythrocytes and Leucocytes are eliminated from plasma)
- Optimal PRP CONCENTRATION



TIME SAVER

- Efficient Processing Time - 2 Minute Single Spin = Time Savings

